

***Baxter***

Baxter Healthcare Ltd., Wallingford Road, Compton, Nr. Newbury,  
Berkshire RG20 7QW Tel. 01635 206000 Fax. 01635 206373

Baxter, FEIBA and BAXJECT are trademarks of Baxter International Inc,  
its subsidiaries or affiliates.

[www.baxterhealthcare.co.uk](http://www.baxterhealthcare.co.uk) [www.haemophiliacare.co.uk](http://www.haemophiliacare.co.uk)

Date of preparation: August 2009. ADV09/2572B



as easy as **1, 2, 3**

A guide to  
reconstituting  
FEIBA therapy

**FEIBA** 

The FEIBA logo consists of the word "FEIBA" in a bold, white, sans-serif font. To the right of the text is a graphic element consisting of three curved lines that sweep upwards and to the right, resembling a stylized swirl or a drop.

## Handy hints

- When withdrawing solution for infusion; the purple side of the BAXJECT II Hi-Flow device will be on top and the word 'Baxter' will be the right way up.
- When withdrawing solution (even if you go wrong); the liquid should always be in the top vial.
- Your syringe is used to provide pressure; either to withdraw solution or to inject air if you go wrong (see the Troubleshooting section for details).